

DESCRIPTION

MEDNA's single-use fluid handling bags are composed of multi-layer HDPE plastic film. They are designed to transfer and store your products. Made from proprietary co-extruded film with ULDPE fluid contacting layer. Assembled using Flexelene™ tubing with 2 male and 2 female luer connectors for easy operator connections. Bags are placed inside a protective packaging. Bags are not sterilized.

SPECIFICATION

10L size

4 ports: 1/4" x 1/4" x 1/4" x 1/4"

Tubing length: 10"

Tubing Diameter:

- O.D: 3/8"
- I.D: 1/4"

2 male luer lock to barb connector

2 female luer lock to barb connector

Supplied with tubing clamps