

DESCRIPTION

The TR-40 Transfection Reagent allows for high-yield protein production in mammalian cell system. In CHO cells, when combined with P54 Feed 1 (Cat # C4611), P54 Feed 2 (Cat # C4612) and P54 Enhancer Solution (Cat # C4613) it provides up to 10-fold improvement to fed batch titers versus batch cultures.

The TR-40 Transfection Reagent allows for high-yield protein production in HEK-293 cells when combined with P40 Feed AC (Cat # C4541) and P40 Feed B (Cat # C4542).

CONDUCTIVITY

16.0 – 19.0 mS/cm

pH

7.0 – 7.3

STERILITY

Negative

STORAGE

Store at 2 - 8°C protected from light.