

DESCRIPTION

P54 Enhancer Solution (20X) is a high protein yield, liquid chemical supplement. It is an animal-origin free partial feed formulation. When combined with P54 Feed 1 (Cat # C4611) and P54 Feed 2 (Cat # C4612) it provides up to 10 fold improvement to fed-batch titers versus unfed batch cultures.

RECOMMENDATION OF FEEDING

One 5% (v/v) addition of Enhancer Solution should be added to transfected cell cultures on the day of transfection.

STERILITY

Negative

STORAGE

Store at 2-8°C protected from light.