

DESCRIPTION

This product is the HRP-Conjugated version of MD10-Monoclonal Anti His-tag IgG (H1307), and it is intended to be used in western blot, ELISA and IHC to detect His-tagged proteins. The conjugates are based on a new method of conjugation to allow high sensitivity and an excellent signal-to-noise ratio.

PRODUCT ADVANTAGES

- High sensitivity
- Excellent signal-to-noise ratio

SUGGESTED WORKING DILUTION RANGE

Chromogenic substrates: 1:250-1:5000

Chemiluminescent substrates: 1:2000-1:100,000

STORAGE

Store at 2-8°C