

DESCRIPTION

Baculovirus particles (BVPs) can be used to measure polyspecificity of therapeutic antibody candidates in in vitro assays. This assesses the viability of an antibody as a therapy early on. Pre-coated BVP ELISA plates are ready to use 96-well strip-well plates coated with BVP and immobilized on the well surface. The plates are supplied in a dry format for easy storage and shipment.

STORAGE

Plates are stable at 20°C for up to twelve months from date of receipt.

Keep plate dry, store unused strips in the packing bag with desiccant.

DATA EXAMPLES

Data shown below is a comparison of plates freshly coated with BVP overnight vs pre-coated BVP dried plates. The ELISA was done using strip-wells on the same plate.

Note: Data for reference purpose only and may differ from actual assays run by end users.

