

DESCRIPTION

Human IgG1 Poly-Specificity Control Antibody is a recombinant antibody recognizing BV particles. It gives poly-specific activity in analytics and molecular applications.

Isotype: IgG1

FORMULATION

100 mM HEPES, 100 mM NaCl, 50 mM NaOAc, pH 6.

STORAGE

Antibody is stable at 2-8°C for up to twelve months from date of receipt.

Keep frozen for prolonged storage.

Data Examples

Data shown below is a typical BVP ELISA using Human IgG1 Poly-Specificity Control Antibody with 6x above BVP score.

Note: Data for reference purpose only and may differ from actual assays run by end users.

