

DESCRIPTION


SE-HPLC size standard is a calibration standard for size exclusion columns used with HPLC in bioanalytics. The mixture includes five proteins represent IgM, BsAb, IgG, Fc and Lysozyme respectively.

The standard can be used with most size exclusion HPLC columns. The standard is supplied as a set of five vials of protein mixture for qualitative analysis only. A 60 uL of standard mixture in each vial is sufficient for 5 injections.

Molecular Weight of Standard Proteins

Proteins MW (kDa)

IgM (peak 1)	900
BsAb (peak 2)	250
IgG (peak 3)	150
Fc (peak 4)	50
Lysozyme (peak 5)	14



STORAGE

Keep frozen.