DESCRIPTION

Human IgG2 Poly-Specificity Negative Control Antibody is a recombinant antibody does not recognize baculovirus particles (BVP). It is used as a negative control in poly-specific activity tests.

Isotype: IgG2

FORMULATION

100 mM HEPES, 100 mM NaCl, 50 mM NaOAc, pH 6.

STORAGE

Antibody is stable at 2-8°C for up to twelve months from date of receipt.
Keep frozen for prolonged storage.

Data Examples

Data shown below is a typical BVP ELISA using Human IgG2 Poly-Specificity Negative Control Antibody with very low affinity to BVP. Positive control antibody (H1801) was probed on the same plate as reference.

Note: Data for reference purpose only and may differ from actual assays run by end users.

